















# **ECG ELECTRODES CATALOGUE**

# ECG ELECTRODES CATALOGUE

REFERENCIA REFERENCE	DESCRIPCIÓN DESCRIPTION	USO USE	PRESENTACIÓN PACKAGING	
99.830.13	Electrodo de corchete 32 x 36 mm. foam con gel sólido adhesivo conductor 32 x 36 mm. foam electrode, press stud connection with adhesive conductor solid gel	ECG en reposo Corta duración Rest ECG Short-term	Caja 600 unidades 12 sobres de 50 unid. 600 units box 12 bags of 50 units	
99.830.26	Electrodo de corchete 32 x 36 mm. foam con gel sólido adhesivo conductor 32 x 36 mm. foam electrode, press stud connection with adhesive conductor solid gel	ECG en reposo Corta duración Rest ECG Short-term	Caja 2000 unidades 20 sobres de 100 unid. 2000 units box 20 bags of 100 units	
99.830.18	Electrodo de corchete 36 x 40 mm. foam con gel sólido adhesivo conductor 36 x 40 mm. foam electrode, press stud connection with adhesive conductor solid gel	ECG en reposo Corta duración Rest ECG Short-term	Caja 600 unidades 20 sobres de 30 unid. 600 units box 20 bags of 30 units	
99.830.11	Electrodo de corchete Ø 50 mm. foam con gel sólido adhesivo conductor Ø 50 mm. foam electrode, press stud connection with adhesive conductor solid gel	ECG en reposo Corta duración Rest ECG Short-term	Caja 600 unidades 20 sobres de 30 unid. 600 units box 20 bags of 30 units	
99.830.12	Electrodo de corchete Ø 30 mm. foam con gel sólido adhesivo conductor Ø 30 mm. foam electrode, press stud connection with adhesive conductor solid gel	Infantil, ECG en reposo Corta duración Rest ECG for children Short-term	Caja 600 unidades 20 sobres de 30 unid. 600 units box 20 bags of 30 units	
99.830.15	Electrodo de tela, corchete de alta sensibilidad Ø 50 mm. gel sólido adhesivo conductor Cloth electrode, press stud connection, high sensitivity Ø 50 mm. adhesive conductor solid gel	Holter Larga monitorización Holter Long-term screening	Caja 600 unidades 20 sobres de 30 unid. 600 units box 20 bags of 30 units	
99.830.16	Electrodo de corchete Clear Tape microperforado Ø 50 mm gel sólido adhesivo conductor. Clear tape microperforated press stud electrode Ø 50 mm. adhesive conductor solid gel.	Holter Larga monitorización Holter Long-term screening	Caja 600 unidades 20 sobres de 30 unid. 600 units box 20 bags of 50 units	

# ECG ELECTRODES CATALOGUE

REFERENCIA REFERENCE	DESCRIPCIÓN DESCRIPTION	USO USE	PRESENTACIÓN PACKAGING	
99.832.01	Electrodo de corchete 57x38mm. foam y gel sólido. Sistema especial fijación latiguillo para evitar artefactos. <i>57x38mm. foam electrode, press-stud connection and solid gel.</i> <i>Special lead fixing system to avoid artifacts.</i>	Prueba de esfuerzo Holter Stress test Holter	Caja 600 unidades 20 sobres de 30 unid. <i>600 units box</i> <i>20 bags of 30 units</i>	
99.832.05	Electrodo de corchete Clear Tape microperforado Ø 50 mm gel sólido adhesivo conductor. Radiotransparente <i>Clear tape microperforated press stud electrode</i> <i>Ø 50 mm. adhesive conductor solid gel. Radiotransparent</i>	Holter Larga monitorización Holter Long-term screening	Caja 600 unidades 20 sobres de 30 unid. <i>600 units box</i> <i>20 bags of 50 units</i>	
99.832.06	Electrodo de corchete Foam 50x48 mm. gel sólido adhesivo conductor. Radiotransparente <i>Foam press stud electrode 50x48 mm.</i> <i>Ø 50 mm. adhesive conductor solid gel. Radiotransparent</i>	ECG en reposo Corta duración Rest ECG Short-term	Caja 600 unidades 20 sobres de 30 unid. <i>600 units box</i> <i>20 bags of 30 units</i>	
99.832.15	Electrodo rectangular 44x32 mm. foam conexión directa banana 4mm. centrada y gel sólido <i>44x32 mm. foam rectangular electrode</i> <i>centered direct 4 mm. banana plug and solid gel</i>	ECG en reposo Corta duración Rest ECG Short-term	Caja 450 unidades 15 sobres de 30 unid. <i>450 units box</i> <i>15 bags of 30 units</i>	
99.832.40	Electrodo rectangular 22 x 34 mm. con lengüeta (tipo TAB) Impregnado con gel sólido adhesivo conductor <i>22 x 34 rectangular electrode with a reed (TAB type)</i> <i>Impregnated with adhesive conductor solid gel</i>	ECG en reposo Corta duración Rest ECG Short-term	Caja 500 unidades 5 sobres de 100 unid. <i>500 units box</i> <i>5 bags of 100 units</i>	



**ENQUIRY: +65 6872 3660**

**EMAIL: [INFO@SCANMED.COM.SG](mailto:INFO@SCANMED.COM.SG)**

**WEBSITE: [SCANMED.COM.SG](http://SCANMED.COM.SG)**